

Procedural Applications

General Surgery

Inguinal Hernia Repair (SS110-30/40/50)
Thyroidectomy (SS110-30/50)
Appendectomy (SS110-70/90)
Splenectomy (SS110-110/140/170)
Distal Pancreatectomy (SS110-110/140/170)
Pancreaticoduodenectomy (SS110-210/260)

Breast Surgery

Sentinel Lymph Node Biopsy (SS110-20/30/40/50)
Lumpectomy (SS110-30/40/50)
Mastectomy (SS110-50/70/90)

Colon and Rectal Surgery

Laparoscopic Colorectal Cancer Resection (SS110-50)
Lap Colectomy (SS110-110/140/170/210/260/310)

Gynecology & Obstetrics

Postpartum Tubal Ligation (SS110-20/30/40)
Bilateral Salpingo-oophorectomy (SS110-30/40/50)
Hysterectomy (SS110-50)
Mini-laparotomy (SS110-50/70/90)
Myomectomy (SS110-50/70/90)
Total Abdominal Hysterectomy (SS110-50/70/90/110/140)
Cesarean Section (SS110-110/140/170/210/260/310)

Bariatric Surgery

Sleeve Gastrectomy (SS110-20/30)
Laparoscopic Gastric Bypass (SS110-40/50)
Open gastric bypass (SS110-110/140/170)

Cardiothoracic Surgery

Video-assisted Thoracoscopic Surgery (SS110-20/30/40/50)
Mitral Valve Repair or Replacement (SS110-50/70/90)
Thoracotomy (SS110-50/70/90)

WOUND PROTECTOR



- 360-degree Protection
- 360-degree Atraumatic Retraction

S

Small Size



Specification	Product Dimensions			Incision size (cm)
	Outer Ring Diameter	Internal ring diameter	Height	
SS110-30	45	40	70	1~2.5
SS110-40	55	50	80	1.5~3
SS110-50	60	70	90	2~4



360° Comprehensive Incision Protection

- Reduces the risk of incision infection.
- Helps prevent tumor seeding and cesarean section endometriosis.
- Maintains a moist wound environment to promote healing.

M

Medium Size



Specification	Product Dimensions			Incision size (cm)
	Outer Ring Diameter	Internal Ring Diameter	Height	
SS110-70	80	90	120	2.5-6
SS110-90	100	110	180	3-7
SS110-110	120	130	150	5-9
SS110-140	150	160	180	7-12

- Easy Operation with a Rigid Retraction Ring
- Suitable for Patients with Different BMI Index
- Accommodates a Wide Range of Specialties
- Competence in Pulling Organs out of the Body for Surgery



- Prevention of Secondary Endometriosis after C-Section
- Shorter Recovery Time
- Smaller Postoperative Scar
- Perfect artificial Birth Canal

L

Large Size



Specification	Product Dimensions			Incision size (cm)
	Outer Ring Diameter	Internal Ring Diameter	Height	
SS110-170	180	190	220	9-140
SS110-210	220	230	280	11-17
SS110-260	270	280	320	22-27

- Adjustable Protective Sleeve
- Minimal Incision Size to 1cm
- Rapid and Effortless Setup

